

CAUTION

It is important that you assist residents in choosing the Personal Help Button (PHB) configuration most suitable for them. Please consider the following:

Pendant Styles

The Philips PHB can be worn around the neck as a pendant by using a neck cord. This helps make the button accessible to both of the senior's hands. **Caution: The neck cord is not designed to break away. Therefore, it can pose a choking risk, including the possibility of death and serious injuries. This may apply to wearers in wheelchairs, using walkers, using beds with guard rails, or who might encounter other protruding objects upon which the cord can become tangled. Wearers for whom this is a concern may wish to consider the Wrist style button.**

Wrist Style

The Philips PHB can also be worn around the wrist with a wristband. **Caution: In circumstances where the user's arm that is not wearing the PHB becomes immobilized (examples: due to stroke or arm being trapped under a heavy object), a risk is that the user will be unable to press the button.**

Please consult your Philips representative regarding which Personal Help Button configuration is most appropriate for your residents.